


Proloquo2Go (PLQ2Go)

1) How are words organized on Proloquo2Go?

Words are color coded (orange outline/background - pronouns, pink outline/background - verbs, green outline/background - prepositions, blue outline/background - adjectives, yellow outline/background - nouns). Many words are organized into “folders”. The folders are mostly category based (e.g. clothes, food). The button locations are organized by some elements of motor planning, which means words stay in approximately the same places on the screen on most pages (e.g. pronouns on the left side). The symbols on PLQ are called SymbolStix. Grammar features, such as verb tenses, are accessed by holding down a verb button.

2) How do I search for a word I can't find on Proloquo2Go?

Click the Tools icon on the bottom of the screen



Then click the magnifying glass.



Type the word you want to find and you will see the “pathway” to find it.

3) What if I want to add words to the system (e.g. sibling names, pet names, specific places we visit, words about specific activities)?

Contact your student’s Speech-Language Pathologist to plan how and where to add the words in the system.

4) How do I access the instruction manual for Proloquo2Go?

<https://www.assistiveware.com/support/proloquo2go/overview/>



5) Is there a training video for Proloquo2Go?

Practical Tips for Setting up PLQ2Go (15 minutes)

<https://www.youtube.com/watch?v=5X3qL8eH21o>



6) What does modeling on Proloquo look like?

Modeling with AAC (2 min, 27 seconds)

<https://www.youtube.com/watch?v=LqeAipcciVQ>



Screen map of Proloquo2Go

