

TouchChat with Word Power (WP)

1) How are words organized on Touchchat with WordPower?

Words are organized into “folders”. Folders are organized grammatically, into categories. The words are also predictive (logical next words and some word completion). For example, verbs automatically change to make grammatically correct sentences (e.g. “drink” turns to “drinking” after “am” is selected so the student automatically says “am” + “drinking” instead of “am” + “drink”). Verb ending (e.g. -ed, -ing) also show up on the next page after a verb is selected. The symbols on TouchChat are called SymbolStix.

2) How do I search for a word I can't find on TouchChat with WP?

Go to “ABC/123” button and then “Find Word” button

How to Use Word Finder (45 seconds)

Video: <https://www.youtube.com/watch?v=H3iO1NOHPPs>



3) What if I want to add words to the system (e.g. sibling names, pet names, specific places we visit, words about specific activities)?

Contact your student's Speech-Language Pathologist to plan how and where to add the words in the system.

4) How do I access the instruction manual for TouchChat Word Power?

<https://touchchatapp.com/support?page=1>



5) Instruction video for TouchChat

TouchChat WordPower60 Basic A Quick Tour (8 min)

<https://www.youtube.com/watch?v=zysnAewI2to>



6) What does modeling look like on TouchChat?

AAC Modeling I Know A Lot book with TouchChat (7 mins)

<https://www.youtube.com/watch?v=thk6xj55Q0c>

